

## CD274

### Research Grade Socazolimab

<b>Catalog No.</b>	DHJ70117A DHJ70117B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	STI-1014, STI-A1014, ZKAB-001, ZKAB001 CAS: 2305043-30-3		
<b>Description:</b>	Research Grade Socazolimab is a humanized monoclonal antibody that targets human Programmed cell death 1 ligand 1 (CD274, B7-H1) Socazolimab is licensed from Sorrento for the Greater China Territory by Lee's Pharma.		
<b>UniProt ID (target):</b>	Q9NZQ7		
<b>Specificity:</b>	Human PD-L1 (CD274, B7-H1)		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	XtenCHO cell line		
<b>Isotype:</b>	Human IgG1 lambda 2		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Applications:</b>	Functional studies		
<b>INN:</b>	Socazolimab		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

