

ANGPTL3

Research Grade Evinacumab

Catalog No.	DHJ88801A DHJ88801B	Quantity:	100 µg 1.0 mg
Alternate Names:	REGN 1500, CAS: 1446419-85-7		
Description:	Research Grade Evinacumab is a human monoclonal antibody directed against human Angiotensin-related protein 3 (Angiotensin-5)		
UniProt ID (target):	Q9Y5C1		
Target:	Human Angiotensin-related protein 3 (Angiotensin-5)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG4 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Evinacumab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

