

## FCRL5

### Research Grade Anti-Human CD307e/FCRL5 (RG7598)

<b>Catalog No.</b>	DHJ94702A DHJ94702B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Fc receptor-like protein 5, Fc receptor homolog 5, FcRH5, RG-7598, CD307e, Immune receptor translocation-associated protein 2, BXMAS1		
<b>Description:</b>	Research Grade Anti-Human CD307e/FCRL5 (RG7598) is a human monoclonal antibody expressed in mammalian cells which targets the tumor-associated antigen (TAA) human Fc receptor-like protein 5		
<b>UniProt ID (target):</b>	Q96RD9		
<b>Target:</b>	Human Fc receptor-like protein 5		
<b>Concentration:</b>	~1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Human IgG1		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.001 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Clone:</b>	RG7598		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or aliquot the stock solution and store at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences®**  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)