

FZD10

Research Grade Tabituximab

| | | | |
|---------------------------------|---|------------------|------------------|
| Catalog No. | DHJ94801A DHJ94801B | Quantity: | 100 µg 1.0 mg |
| Alternate Names: | ONC-005, OTSA101, OTSA101-DTPA-90Y, OTSA101-DTPA, OTSA101-SS01, CAS: 1644134-10-0 | | |
| Description: | Research Grade Tabituximab is a chimeric monoclonal antibody which targets human Frizzled-10 (CD350) | | |
| UniProt ID (target): | Q9ULW2 | | |
| Target: | Human Frizzled-10 (CD350) | | |
| Concentration: | 1.0 mg/ml, lot specific | | |
| Source: | XtenCHO | | |
| Isotype: | Chimeric IgG1 kappa | | |
| Formulation: | 0.01M PBS buffer, pH 7.4 | | |
| Purity: | > 95% by SDS-PAGE | | |
| Endotoxin Level: | ≤ 0.01 EU/µg by LAL analysis | | |
| Purification: | Protein A/G affinity chromatography | | |
| INN: | Tabituximab | | |
| Applications: | Functional studies | | |
| Storage & Stability: | Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles. | | |

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

