

Research Grade Satumomab

Catalog No.	DHK10001A DHK10001B	Quantity:	100 µg 1.0 mg
Alternate Names:	B72.3, Indium (111In) satumomab pendetide, OncoScint CR/OV, OncoScint CR103, CAS: 144058-40-2		
Description:	Research Grade Satumomab is a murine monoclonal antibody directed against human tumor-associated glycoprotein 72 (TAG-72)		
Target:	Human tumor-associated glycoprotein 72 (TAG-72)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Mouse IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Satumomab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.