

## SIGLEC15

### Research Grade Anti-Human SIGLEC15/CD33L3 (DS1501) mAb

<b>Catalog No.</b>	DHK13704A DHK13704B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Sialic acid-binding Ig-like lectin 15, CD33L3, Siglec-15, SIGLEC15, CD33 antigen-like 3		
<b>Description:</b>	Research Grade Anti-Human SIGLEC15/CD33L3 (DS1501) is a humanized monoclonal antibody which targets human Sialic acid-binding Ig-like lectin 15 (Siglec-15, CD33L3)		
<b>UniProt ID (target):</b>	Q6ZMC9		
<b>Specificity:</b>	Human Siglec-15 (CD33L3)		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	XtenCHO		
<b>Isotype:</b>	Human IgG1 kappa (silent)		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A or G affinity chromatography		
<b>Applications:</b>	Functional studies		
<b>Clone:</b>	DS1501		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

