

## TTR

### Research Grade Coramitug

<b>Catalog No.</b>	DHK25501A DHK25501B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Transthyretin, Prealbumin, TBPA, TTR, PALB, ATTR, Transthyretin amyloidosis, CAS: 2626841-92-5		
<b>Description:</b>	Research Grade Coramitug is a chimeric monoclonal antibody affinity purified from cell culture supernatant. The investigational humanized monoclonal antibody coramitug was designed to target and clear the non-native transthyretin aggregates (misTTR depleter) associated with disease pathology that underlies both wtATTR and hATTR, without affecting the native, or normal tetrameric form of the protein (TTR).		
<b>UniProt ID (target):</b>	P02766-A		
<b>Specificity:</b>	Human Transthyretin		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Chimeric humanized IgG1 kappa		
<b>Formulation:</b>	Sterile filtered 0.01M PBS, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A or G affinity chromatography from cell culture supernatant		
<b>INN:</b>	Coramitug		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)