

CRLF2

Research Grade Verekitug

Catalog No.	DHK27201A DHK27201B	Quantity:	100 µg 1.0 mg
Alternate Names:	Cytokine receptor-like factor 2, IL-XR, Thymic stromal lymphopoietin protein receptor, TSLPR, CAS: 2768386-15-6		
Description:	Research Grade Verekitug is a human monoclonal antibody directed against human Cytokine receptor-like factor 2		
UniProt ID (target):	Q9HC73		
Target:	Human Cytokine receptor-like factor 2		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Human IgG1 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Verekitug		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

