

CRLF2

Research Grade Anti-Human TSLPR (RG-7258) mAb

Catalog No.	DHK27202A DHK27202B	Quantity:	100 µg 1.0 mg
Alternate Names:	Cytokine receptor-like factor 2, IL-XR, Thymic stromal lymphopoietin protein receptor, TSLPR, RG7258, RG 7258		
Description:	Research Grade Anti-Human TSLPR (RG-7258) is a humanized monoclonal antibody directed against human Cytokine receptor-like factor 2, purified by affinity chromatography.		
UniProt ID (target):	Q9HC73		
Target:	Human Cytokine receptor-like factor 2		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Human IgG1 kappa		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.001 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
INN:	RG-7258		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences®
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com