

HLA-DQB1

Research Grade Humanized DQB1-FL8

Catalog No.	DHM01501A DHM01501B	Quantity:	100 µg 1.0 mg
Alternate Names:	DQB1-FL8, LG11		
Description:	Research Grade Humanized DQB1-FL8 is a monoclonal antibody directed against human HLA class II histocompatibility antigen, DQ beta 1 chain		
UniProt ID (target):	P01920		
Target:	Human HLA class II histocompatibility antigen, DQ beta 1 chain		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1		
Formulation:	Sterile-filtered 0.01 M PBS, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.