

S

Research Grade Adintrevimab

Catalog No.	DVV00315A DVV00315B	Quantity:	100 µg 1.0 mg
Alternate Names:	ADG-2, ADG-20 CAS: 2516243-54-0, SARS-like Coronaviruses		
Description:	Research Grade Adintrevimab is a human monoclonal antibody that targets Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2) Spike glycoprotein		
UniProt ID (target):	P0DTC2-2		
Specificity:	SARS-CoV-2 Spike glycoprotein EC ₅₀ = 37.14 with RBD by indirect ELISA		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO cell line		
Isotype:	Human IgG1 lambda2		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies		
INN:	Adintrevimab		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com