

F

Research Grade HNK20

Catalog No.	DVV02805A DVV02805B	Quantity:	100 µg 1.0 mg
Alternate Names:	Human Respiratory Syncytial Virus (HRSV)		
Description:	Research Grade HNK20 is a mouse monoclonal antibody directed against Fusion glycoprotein F0 of Human respiratory syncytial virus A (HRSV-A).		
UniProt ID (target):	P03420		
Target:	Fusion glycoprotein F0 of HRSV-A		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Murine IgGA		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

