

S

Research Grade Tuvirumab

Catalog No.	DVV03001A DVV03001B	Quantity:	100 µg 1.0 mg
Alternate Names:	OST 577, CAS: 138660-97-6, Hepatitis B Virus (HBV-A)		
Description:	Research Grade Tuvirumab is a human monoclonal antibody directed against Large envelope protein of Hepatitis B Virus (HBV-A)		
UniProt ID (target):	P03142		
Target:	Large envelope protein of Hepatitis B Virus (HBV-A)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 lambda		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Tuvirumab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.