

## stx2 A-subunit Research Grade Setoxaximab

<b>Catalog No.</b>	DXX00701A DXX00701B	<b>Quantity:</b>	100 µg 1.0 mg
<b>Alternate Names:</b>	Shigamabs (drug name for the combination of pritoxaximab and setoxaximab), caStx2, CAS: 1351470-17-1		
<b>Description:</b>	Research Grade Setoxaximab is a chimeric monoclonal antibody directed against Shiga toxin 2 A subunit (rRNA N-glycosylase) of <i>E. coli</i> .		
<b>UniProt ID (target):</b>	Q8XBV2		
<b>Target:</b>	Shiga toxin 2 A subunit (rRNA N-glycosylase) of <i>E. coli</i> .		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	XtenCHO		
<b>Isotype:</b>	Chimeric IgG1 kappa		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Endotoxin Level:</b>	≤ 0.01 EU/µg by LAL analysis		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>INN:</b>	Setoxaximab		
<b>Applications:</b>	Functional studies		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

