

Research Grade Edobacomab

Catalog No.	DXX02202A DXX02202B	Quantity:	100 µg 1.0 mg
Alternate Names:	E5, XMMEN-0E5, CAS: 141410-98-2, Bacteria		
Description:	Research Grade Edobacomab is a mouse monoclonal antibody directed against the Lipid A domain of endotoxin in Gram-negative bacteria. Edobacomab is under investigation for the treatment of sepsis caused by Gram-negative bacteria.		
Target:	Gram-negative bacterial Endotoxin (Lipid A domain)		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Mouse IgM		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Edobacomab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.