

HSP90

Research Grade Efungumab

Catalog No.	DXX02301A DXX02301B	Quantity:	100 µg 1.0 mg
Alternate Names:	HSP90mab, Mycograb, CAS: 762260-74-2		
Description:	Research Grade Efungumab is a human monoclonal antibody directed against Heat shock protein 90 homolog of <i>Candida albicans</i> .		
UniProt ID (target):	P46598		
Target:	Heat shock protein 90 homolog of <i>Candida albicans</i> .		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human scFv-heavy-kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Efungumab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

