

pcrV & psl

Research Grade Gremubamab

Catalog No.	DXX03401A DXX03401B	Quantity:	100 µg 1.0 mg
Alternate Names:	Bispecific, MEDI3902, CAS: 1800381-36-5		
Description:	Research Grade Gremubamab is a bispecific humanized monoclonal antibody directed against Type III secretion protein PcrV and Polysaccharide synthesis locus (Psl) of <i>Pseudomonas aeruginosa</i> , designed for the treatment of ventilator-associated pneumonia (VAP).		
UniProt ID (target):	O30527		
Target:	Type III secretion protein PcrV and Polysaccharide synthesis locus (Psl) of <i>Pseudomonas aeruginosa</i>		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human VH-CH1-VH-V-KAPPA-CH2-CH3 (IgG1 kappa)		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/µg by LAL analysis		
Purification:	Protein A/G affinity chromatography		
INN:	Gremubamab		
Applications:	Functional studies		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com