

IGHE

Recombinant Rat IgE Protein

Catalog No.	ERJ92701A ERJ92701B	Quantity:	100 µg 1.0 mg
Alternate Names:	Immunoglobulin heavy constant epsilon, Ig epsilon chain C region		
Description:	Full length recombinant rat IgE was expressed in mammalian cells		
UniProt ID:	P01855		
Species:	<i>Rattus norvegicus</i>		
Source:	Mammalian cells		
Molecular Weight:	190 kDa predicted		
Concentration:	1.0 mg/ml, lot specific		
Formulation:	PBS pH 6.0 containing 0.2M Arg		
Purity:	> 90% as determined by SDS-PAGE		
Storage & Stability:	Store as supplied at -20°C to -80°C for up to 1 year. Upon reconstitution, prepare working aliquots and store at -20°C to -80°C for up to 3 months. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

