

CU43 Protease (Human IgG Specific)

Catalog No.	EXX25001A EXX25001B	Quantity:	100 µg 1.0 mg
Alternate Names:	CU43 Protease, CU43		
Description:	Recombinant CU43 Protease is cloned from <i>Corynebacterium ulcerans</i> and expressed in mammalian cells (Note: YXX25001 expressed in <i>E. coli</i>)		
Species:	<i>Corynebacterium ulcerans</i>		
Source:	Expressed in Mammalian cells		
Formulation:	Sterile-filtered PBS, pH 7.4		
Molecular Weight:	70.94 kDa		
Purity:	> 90 % as determined by SDS-PAGE.		
Applications:	Proteomics		
Storage & Stability:	Stable for up to 1 year from date of receipt at -20°C to -80°C or store working aliquots at -20°C to -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

