

BTN3A1

Mouse Anti-Human CD277/BTN3A1 (103.2#) mAb, PE

Catalog No.	FHA13512A FHA13512B	Quantity:	50 tests 100 tests
Alternate Names:	Butyrophilin subfamily 3 member A1, BTN3A1, BTF5, CD277, hu103.2		
Description:	Mouse monoclonal Anti-Human CD277/BTN3A1 (103.2#) was purified by affinity chromatography and conjugated to phycoerythrin (PE).		
UniProt ID (target):	O00481		
Target:	Human Butyrophilin subfamily 3 member A1 (CD277)		
Concentration:	0.5 mg/ml		
Source:	Mammalian cell line		
Isotype:	Mouse IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4, 0.2% BSA, 0.05% Proclin 300.		
Purity:	> 95% by SDS-PAGE		
Clone:	103.2#		
Purification:	Protein A/G affinity chromatography		
Applications:	Flow cytometry		
Conjugate:	Phycoerythrin		
Storage & Stability:	Store at 2-8°C for 1 year. Protect from light. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

