

GDF11

Mouse Anti-Human GDF11/BMP11 (SAA0431) mAb

Catalog No.	FHB63020A FHB63020B	Quantity:	100 µg 1.0 mg
Alternate Names:	Growth/differentiation factor 11, GDF-11, Bone morphogenetic protein 11, BMP-11		
Description:	Mouse Anti-Human GDF11/BMP11 (SAA0431) mAb is purified from cell culture supernatant.		
UniProt ID (target):	O95390		
Target:	Human GDF-11		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4, 0.09% Sodium azide.		
Purity:	> 95% by SDS-PAGE and SEC-HPLC		
Clone:	SAA0431		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	ELISA, FCM		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

