

CD4

Mouse Anti-Human CD4 (Clone OKT4A) mAb

Catalog No.	FHB95930A FHB95930B	Quantity:	50 µg 100 µg
Alternate Names:	T-cell surface glycoprotein CD4, T-cell surface antigen T4/Leu-3, CD4		
Description:	Mouse Anti-human CD4 (Clone OKT4A) was affinity purified from cell culture supernatant and targets Human T-cell surface glycoprotein CD4		
UniProt ID (target):	P01730		
Target:	Human T-cell surface glycoprotein CD4		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Applications:	Functional studies, FCM, IHC, IF		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

