

CD4

Mouse Anti-Human CD4 (Clone OKT4A) mAb, PE

Catalog No.	FHB95932A FHB95932B	Quantity:	50 tests 100 tests
Alternate Names:	T-cell surface glycoprotein CD4, T-cell surface antigen T4/Leu-3, CD4		
Description:	Mouse Anti-human CD4 (Clone OKT4A) was affinity purified from cell culture supernatant, targets Human T-cell surface glycoprotein CD4., and conjugated to phycoerythrin.		
UniProt ID (target):	P01730		
Target:	Human T-cell surface glycoprotein CD4		
Concentration:	0.5 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.		
Purity:	> 95% by SDS-PAGE		
Clone:	OKT4A		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM, IF		
Conjugate:	Phycoerythrin (PE)		
Storage & Stability:	Store at 2-8°C for up to 1 year. Protect from light. Do not freeze		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

