

ITGB3

Mouse Anti-Human CD61/ITGB3 (Clone SZ-21) mAb

Catalog No.	FHC13020A FHC13020B	Quantity:	50 µg 100 µg
Alternate Names:	Integrin beta-3, ITGB3, Platelet membrane glycoprotein IIIa, GPIIIa, CD61		
Description:	Mouse Anti-human CD61/ITGB3 (Clone SZ-21) was affinity purified from cell culture supernatant.		
UniProt ID (target):	P05106		
Target:	Human Integrin beta-3		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG1, kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	SZ-21		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	FCM		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.