

KLK3

Anti-Human PSA/KLK3 (SAA0506) mAb

Catalog No.	FHC24210A FHC24210B	Quantity:	50 µg 100 µg
Alternate Names:	Prostate-specific antigen, PSA, Gamma-seminoprotein, Semin, Kallikrein-3, P-30 antigen, Semenogelase		
Description:	Anti-Human PSA/KLK3 (SAA0506) monoclonal antibody is purified from cell culture supernatant by affinity chromatography.		
UniProt ID (target):	P07288		
Target:	Human prostate-specific antigen		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Human IgG1		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	SAA0506		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	Blocking, ELISA, Flow cytometry		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

