

## GP1BA

### Mouse Anti-Human CD42b/GP1BA (Clone 5G6) mAb

<b>Catalog No.</b>	FHC25710A FHC25710B	<b>Quantity:</b>	50 µg 100 µg
<b>Alternate Names:</b>	Platelet glycoprotein Ib alpha chain, GP-Ib alpha; GPIb-alpha; GPIbA; Glycoprotein Ibalpha, GP-Ib alpha; GPIb-alpha; GPIbA; Glycoprotein Ibalpha, Glycocalicin, CD42b		
<b>Description:</b>	Mouse Anti-Human CD42b/GP1BA (Clone 5G6) mAb is purified from cell culture supernatant by affinity chromatography.		
<b>UniProt ID (target):</b>	P07359		
<b>Target:</b>	Human CD42b		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Mouse IgG1 kappa		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Clone:</b>	5G6		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Applications:</b>	FCM, Western blot		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)