

CD3E

Mouse Anti-Human CD3 (Clone SP34) mAb

Catalog No.	FHC27750A FHC27750B	Quantity:	50 µg 100 µg
Alternate Names:	T-cell surface glycoprotein CD3 epsilon chain, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e		
Description:	Mouse Anti-Human CD3 (Clone SP34) mAb is purified from cell culture supernatant by affinity chromatography.		
UniProt ID (target):	P07766		
Target:	Human CD3		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4 containing 0.09% sodium azide		
Purity:	> 95% by SDS-PAGE		
Clone:	SP34		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

