

ITGA2B & ITGB3

Anti-Human CD41/ITGA2B & CD61/ITGB3 (7E3#) mAb, APC

Catalog No.	FHC33310-APC	Quantity:	100 tests
Alternate Names:	Integrin alpha-IIb, CD41, GPalpha IIb, GPIIb, Platelet membrane glycoprotein IIb, Integrin beta-3, CD61, Platelet membrane glycoprotein IIIa, GPIIIa, Allophycocyanin, APC		
Description:	Anti-Human CD41/ITGA2B & CD61/ITGB3 (7E3#) is a human monoclonal antibody that targets human Integrin alpha-IIb (CD41) and Integrin beta-3 (CD61), conjugated to Allophycocyanin (APC)		
UniProt ID (target):	P08514 & P05106		
Target:	Human Integrin alpha-IIb (CD41) and Integrin beta-3 (CD61)		
Source:	Mammalian cells		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	7E3#		
Purification:	Protein A/G affinity chromatography		
Conjugate:	Allophycocyanin (APC)		
Applications:	FCM		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com