

PTPRC

Mouse Anti-Human CD45/PTPRC (BC8) mAb, APC

Catalog No.	FHC34013A FHC34013B	Quantity:	50 µg 100 µg
Alternate Names:	Receptor-type tyrosine-protein phosphatase C, CD45, T200, Leukocyte common antigen, L-CA		
Description:	Mouse Anti-Human CD45/PTPRC (BC8) mAb, APC is a monoclonal antibody purified from cell culture supernatant by affinity chromatography and conjugated to allophycocyanin for use in flow cytometry.		
UniProt ID (target):	P08575		
Target:	Human Receptor-type tyrosine-protein phosphatase C		
Concentration:	0.5 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG1 kappa		
Formulation:	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
Purity:	> 95% by SDS-PAGE		
Clone:	BC8		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	FCM		
Conjugate:	Allophycocyanin		
Storage & Stability:	Store at 2-8°C for 12 months. Protect from light. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

