

ITGAM

Mouse Anti-Human CD11b/ITGAM (mAb107) mAb, APC

Catalog No.	FHC86613A FHC86613B	Quantity:	50 tests 100 tests
Alternate Names:	Integrin alpha-M, CD11b, CD11 antigen-like family member B, CR-3 alpha chain, Cell surface glycoprotein MAC-1 subunit alpha, Leukocyte adhesion receptor MO1, Neutrophil adherence receptor		
Description:	Mouse Anti-Human CD11b/ITGAM (mAb107) mAb is purified from cell culture supernatant by affinity chromatography and conjugated to allophycocyanin (APC).		
UniProt ID (target):	P11215		
Target:	Human CD11b/ITGAM		
Conjugate:	Allophycocyanin (APC)		
Source:	Mammalian cells		
Isotype:	Mouse IgG1 kappa		
Formulation:	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300.		
Purity:	> 95% by SDS-PAGE		
Clone:	mAb107		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM, Immunoblot		
Concentration:	0.5 mg/ml, lot specific		
Storage & Stability:	Store at 2-8°C for 1 year. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

