

FOLR1

Anti-Human FOLR1 (SAA0118) mAb, PE

Catalog No.	FHD10422A FHD10422B	Quantity:	50 tests 100 tests
Alternate Names:	Folate receptor alpha, FR-alpha, Adult folate-binding protein, FBP, Folate receptor 1, Folate receptor, adult, KB cells FBP, Ovarian tumor-associated antigen MOv18		
Description:	Anti-Human FOLR1 (SAA0118) mAb is a affinity purified humanized monoclonal antibody conjugated to phycoerythrin which targets human Folate receptor alpha (FR-alpha)		
UniProt ID (target):	P15328		
Specificity:	Human Folate receptor alpha (FR-alpha)		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
Purity:	> 95% by SDS-PAGE		
Endotoxin Level:	≤ 0.01 EU/μg by LAL analysis		
Purification:	Protein A or G affinity chromatography		
Clone:	SAA0118		
Conjugate:	Phycoerythrin (PE)		
Applications:	Flow cytometry		
Storage & Stability:	Store at 2-8°C for up to 1 year. Protect from light. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.