

## VEGFA

### Anti-Human VEGFA (SAA0446) mAb, APC

<b>Catalog No.</b>	FHD12613A FHD12613B	<b>Quantity:</b>	50 tests 100 tests
<b>Alternate Names:</b>	Vascular endothelial growth factor A, VEGFA, VEGF, Vascular permeability factor, VEGF-A		
<b>Description:</b>	Anti-Human VEGFA (SAA0446) mAb is a affinity purified mouse monoclonal antibody conjugated to allophycocyanin.		
<b>UniProt ID (target):</b>	P15692		
<b>Specificity:</b>	Human Vascular endothelial growth factor A		
<b>Concentration:</b>	0.5 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Purification:</b>	Protein A or G affinity chromatography from cell culture supernatant		
<b>Clone:</b>	SAA0446		
<b>Conjugate:</b>	Allophycocyanin (APC)		
<b>Applications:</b>	Flow cytometry, ELISA		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year. <b>Protect from light. Do not freeze.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

