

PRLR

Humanized Anti-Human PRLR (Clone h16f) mAb

Catalog No.	FHD18010A FHD18010B	Quantity:	100 µg 1.0 mg
Alternate Names:	Prolactin receptor, PRLR, PR-1594804, h16f,		
Description:	Humanized monoclonal antibody derived from the mouse monoclonal antibody (clone h16f) developed to target the human prolactin receptor. This antibody is the basis for Rolinsatamab (PR-1594804) and antibody drug conjugate (ADC) ABBV-176		
UniProt ID (target):	P16471		
Target:	Human prolactin receptor		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Humanized IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.