

ACVR2A & ACVR2B

Anti-Human ACVR2A & ACVR2B (SAA0125) mAb, PE

Catalog No.	FHD72812A FHD72812B	Quantity:	50 tests 100 tests
Alternate Names:	Activin receptor type-2A, Activin receptor type IIA, ACTR-IIA, ACTRIIA, Activin receptor type-2B, Activin receptor type IIB, ACTR-IIB		
Description:	Anti-Human Activin receptor type-2A & 2B (SAA0125) mAb, PE is a affinity purified humanized monoclonal antibody conjugated to phycoerythrin.		
UniProt ID (target):	P27037 & Q13705		
Specificity:	Human Activin receptor type-2A & 2B		
Concentration:	0.5 mg/ml		
Source:	Mammalian cells		
Isotype:	Human IgG1 lambda		
Formulation:	0.01M PBS buffer, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Clone:	SAA0125		
Conjugate:	Phycoerythrin (PE)		
Applications:	Flow cytometry		
Storage & Stability:	Store at 2-8°C for up to 1 year. Protect from light. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

