

CD38

Anti-Human CD38/ADPRC1 (HuMax-003) mAb, FITC

Catalog No.	FHD80851A FHD80851B	Quantity:	50 tests 100 tests
Alternate Names:	ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1, ADP-ribosyl cyclase 1, ADPRC 1, CD38		
Description:	Anti-Human CD38/ADPRC1 (HuMax-003) is a affinity purified mouse monoclonal antibody conjugated to FITC which targets human ADP-ribosyl cyclase 1/CD38		
UniProt ID (target):	P28907		
Specificity:	Human ADP-ribosyl cyclase 1/CD38		
Concentration:	0.5 mg/ml		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A or G affinity chromatography		
Clone:	HuMax-003		
Conjugate:	FITC		
Applications:	Flow cytometry, ELISA		
Storage & Stability:	Store at 2-8°C for up to 1 year. Protect from light. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

