

## CD38

### Anti-Human CD38/ADPRC1 (HuMax-003) mAb, APC

<b>Catalog No.</b>	FHD80853A FHD80853B	<b>Quantity:</b>	50 tests 100 tests
<b>Alternate Names:</b>	ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1, ADP-ribosyl cyclase 1, ADPRC 1, CD38		
<b>Description:</b>	Anti-Human CD38/ADPRC1 (HuMax-003) is a affinity purified mouse monoclonal antibody conjugated to APC which targets human ADP-ribosyl cyclase 1/CD38		
<b>UniProt ID (target):</b>	P28907		
<b>Specificity:</b>	Human ADP-ribosyl cyclase 1/CD38		
<b>Concentration:</b>	0.5 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Clone:</b>	HuMax-003		
<b>Purification:</b>	Protein A or G affinity chromatography		
<b>Conjugate:</b>	Allophycocyanin (APC)		
<b>Applications:</b>	Flow cytometry, ELISA		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 year. <b>Protect from light. Do not freeze.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

