

FCGR2B

Anti-Human CD32b/FCGR2B (SAA0547) mAb, PE

Catalog No.	FHD97912A FHD97912B	Quantity:	50 tests 100 tests
Alternate Names:	Low affinity immunoglobulin gamma Fc region receptor II-b, IgG Fc receptor II-b, CDw32, Fc-gamma RII-b, Fc-gamma-RIIb, FcRII-b		
Description:	Anti-Human CD32b/FCGR2B (SAA0547) mAb, PE is affinity purified and labeled with phycoerythrin.		
UniProt ID (target):	P31994		
Target:	Human FCGR2B		
Concentration:	0.5 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Human IgG1, kappa		
Formulation:	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Clone:	SAA0547		
Applications:	FCM		
Conjugate:	Phycoerythrin (PE)		
Storage & Stability:	Store at 2-8°C for up to 1 year. Protect from light. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.