

## FCGR2B

### Mouse Anti-Human CD32b/FCGR2B (6G08) mAb

<b>Catalog No.</b>	FHD97920A FHD97920B	<b>Quantity:</b>	50 µg 100 µg
<b>Alternate Names:</b>	Low affinity immunoglobulin gamma Fc region receptor II-b, IgG Fc receptor II-b, CDw32, Fc-gamma RII-b, Fc-gamma-RIIb, FcRII-b		
<b>Description:</b>	Mouse Anti-Human CD32b/FCGR2B (6G08) mAb is affinity purified from cell culture supernatant		
<b>UniProt ID (target):</b>	P31994		
<b>Target:</b>	Human FCGR2B		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Human IgG1 lambda		
<b>Formulation:</b>	0.01M PBS, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Purification:</b>	Protein A/G affinity chromatography from cell culture supernatant		
<b>Clone:</b>	6G08		
<b>Applications:</b>	FCM		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for immediate use (1 week) or in aliquots at -20°C to -80°C for up to 1 year. <b>Avoid repeated freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



**Cell Sciences**<sup>®</sup>  
65 Parker Street  
Unit 11  
Newburyport, MA 01950

Toll Free: 888-769-1246  
Phone: 978-572-1070  
Fax: 978-992-0298

E-mail: [info@cellsciences.com](mailto:info@cellsciences.com)  
Website: [www.cellsciences.com](http://www.cellsciences.com)