

CD96

Mouse Anti-Human CD96 (BA072) mAb

Catalog No.	FHE25710A FHE25710B	Quantity:	50 µg 100 µg
Alternate Names:	T-cell surface protein tactile, Cell surface antigen CD96, T cell-activated increased late expression protein		
Description:	Affinity purified mouse monoclonal antibody that targets Human T-cell surface protein tactile/CD96 for use in flow cytometry.		
UniProt ID (target):	P40200		
Target:	Human T-cell surface protein tactile/CD96		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	BA072		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	FCM		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

