

PMEL

Mouse Anti-Human PMEL17/ME20M (17A9), PE

Catalog No.	FHE28712A FHE28712B	Quantity:	50 tests 100 tests
Alternate Names:	Melanocyte protein PMEL, ME20-M, ME20M, Melanocyte protein Pmel 17, Melanocytes lineage-specific antigen GP100, Melanoma-associated ME20 antigen, Premelanosome protein, Silver locus protein homolog		
Description:	Mouse Anti-Human PMEL17/ME20M (17A9) also reacts with mouse, rat and monkey Melanocyte protein PMEL.		
UniProt ID (target):	P40967		
Target:	Human Melanocyte protein PMEL		
Concentration:	0.5 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
Purity:	> 95% by SDS-PAGE		
Clone:	17A9		
Purification:	Protein A/G affinity chromatography		
Applications:	ELISA, FCM (0.5 µg/test)		
Cross Reactivity:	Mouse, Rat, Monkey		
Conjugate:	Phycoerythrin (PE)		
Storage & Stability:	Store at 2-8°C for 1 year. Protect from light. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

