

ROR1

Human Anti-Human ROR1 (SAA0140) mAb

Catalog No.	FHF90710A FHF90710B FHF90710C	Quantity:	50 µg 100 µg 1.0 mg
Alternate Names:	Inactive tyrosine-protein kinase transmembrane receptor ROR1, mROR1, Neurotrophic tyrosine kinase, receptor-related 1		
Description:	Receptor tyrosine kinase ROR1 plays an essential role in embryogenesis and is overexpressed in many types of malignant tumors. Studies have demonstrated that in oncogenesis ROR1 activates cell survival signaling events, particularly the non-canonical WNT signaling pathway. Antibody-based immunotherapies targeting ROR1 show promise in certain malignancies.		
UniProt ID (target):	Q01973		
Target:	Human ROR1		
Concentration:	1.0 mg/ml, lot specific		
Source:	XtenCHO		
Isotype:	Human IgG1 kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Clone:	SAA0140		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

