

MSLN

Mouse Anti-Human MSLN/Mesothelin (SS1) mAb, FITC

Catalog No.	FHG52041A FHG52041B	Quantity:	50 Tests 100 Tests
Alternate Names:	Mesothelin, CAK1 antigen, Pre-pro-megakaryocyte-potentiating factor, MPF, MSLN		
Description:	Mouse Anti-human Mesothelin (m912) was affinity purified from cell culture supernatant.		
UniProt ID (target):	Q13421		
Target:	Human Mesothelin		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
Purity:	> 95% by SDS-PAGE		
Conjugate:	Fluorescein isothiocyanate (FITC)		
Purification:	Protein A/G affinity chromatography		
Applications:	ELISA, FCM		
Storage & Stability:	Store at 2-8°C for 1 year. Protect from light. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

