

## IL11RA

### Mouse Anti-Human IL11RA (8E4#) mAb

<b>Catalog No.</b>	FHG84920A FHG84920B	<b>Quantity:</b>	50 µg 100 µg
<b>Alternate Names:</b>	Interleukin-11 receptor subunit alpha, L-11R subunit alpha, IL-11R-alpha, IL-11RA		
<b>Description:</b>	Mouse Anti-human IL11RA (8E4#) was affinity purified from cell culture supernatant.		
<b>UniProt ID (target):</b>	Q14626		
<b>Target:</b>	Human Mesothelin		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Mouse IgG2a kappa		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Clone:</b>	8E4#		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Applications:</b>	Blocking, ELISA, FCM		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**