

MICA

Anti-Human MICA (1D5) mAb

Catalog No.	FHH45230A FHH45230B	Quantity:	50 µg 100 µg
Alternate Names:	MHC class I polypeptide-related sequence A, MIC-A		
Description:	Anti-Human MICA (1D5) is a mouse monoclonal antibody expressed in mammalian cells and affinity purified.		
UniProt ID (target):	Q29983		
Specificity:	Human MHC class I polypeptide-related sequence A		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cell line		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Clone:	1D5		
Purification:	Protein A/G affinity chromatography from cell culture supernatant		
Applications:	FCM, ELISA		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

