

VTCN1

Mouse Anti-Human VTCN1/B7-H4 (6H3) mAb, APC

Catalog No.	FHJ14833A FHJ14833B	Quantity:	50 tests 100 tests
Alternate Names:	V-set domain-containing T-cell activation inhibitor 1, B7 homolog 4, B7-H4, B7h.5, Immune costimulatory protein B7-H4, Protein B7S1, T-cell costimulatory molecule B7x		
Description:	Affinity purified Mouse Anti-Human VTCN1/B7-H4 (6H3) monoclonal antibody targets Human V-set domain-containing T-cell activation inhibitor 1 and is conjugated to APC.		
UniProt ID (target):	Q7Z7D3		
Target:	Human V-set domain-containing T-cell activation inhibitor 1		
Concentration:	0.5 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
Purity:	> 95% by SDS-PAGE		
Clone:	6H3		
Purification:	Protein A/G affinity chromatography		
Conjugate:	APC		
Applications:	FCM (0.5µg/test) ELISA, IF		
Storage & Stability:	Store at 2-8°C for 12 months. Protect from light. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

