

LRRC15

Mouse Anti-Human LRRC15/LIB (SAA0773) mAb

Catalog No.	FHJ28810A FHJ28810B	Quantity:	100 µg 1.0 mg
Alternate Names:	Leucine-rich repeat-containing protein 15, LRRC15, Leucine-rich repeat protein induced by beta-amyloid homolog, hLib		
Description:	Monoclonal Mouse Anti-Human LRRC15/LIB (SAA0773) targets human Leucine-rich repeat-containing protein 15.		
UniProt ID (target):	Q8TF66		
Target:	Human Leucine-rich repeat-containing protein 15 (LRRC15)		
Concentration:	1.0 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4		
Purity:	> 95% by SDS-PAGE		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM, ELISA		
Storage & Stability:	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

