

LRRC15

Mouse Anti-Human LRRC15/LIB (SAA0773) mAb, FITC

Catalog No.	FHJ28811A FHJ28811B	Quantity:	50 tests 100 tests
Alternate Names:	Leucine-rich repeat-containing protein 15, LRRC15, Leucine-rich repeat protein induced by beta-amyloid homolog, hLib		
Description:	Monoclonal Mouse Anti-Human LRRC15/LIB (SAA0773), FITC targets human Leucine-rich repeat-containing protein 15 for use in flow cytometry.		
UniProt ID (target):	Q8TF66		
Target:	Human Leucine-rich repeat-containing protein 15 (LRRC15)		
Concentration:	0.5 mg/ml, lot specific		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS buffer, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
Purity:	> 95% by SDS-PAGE		
Conjugate:	FITC		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM, ELISA		
Clone:	SAA0773		
Storage & Stability:	Store at 2-8°C for 12 months. Protect from light. DO NOT FREEZE.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

