

NECTIN2

Mouse Anti-Human CD112/Nectin-2/PVRL2 (SAA0405) mAb, PerCP

Catalog No.	FHJ32514A FHJ32514B	Quantity:	50 tests 100 tests
Alternate Names:	Nectin-2, Herpes virus entry mediator B, HveB, Nectin cell adhesion molecule 2, Poliovirus receptor-related protein 2, PRR2, CD112		
Description:	Monoclonal Mouse Anti-Human CD112/Nectin (SAA0405) which targets human CD112/Nectin-2 is affinity purified and labeled with Peridinin-Chlorophyll-Protein (PerCP)		
UniProt ID (target):	Q92692		
Target:	Human Nectin-2		
Source:	Mammalian cells		
Isotype:	Mouse IgG2a kappa		
Formulation:	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
Purity:	> 95% by SDS-PAGE		
Clone:	SAA0405		
Purification:	Protein A/G affinity chromatography		
Applications:	FCM, ELISA		
Conjugate:	Peridinin-Chlorophyll-Protein (PerCP)		
Storage & Stability:	Store at 2-8°C for 12 months. Do not freeze.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences[®]
65 Parker Street
Unit 11
Newburyport, MA 01950

Toll Free: 888-769-1246
Phone: 978-572-1070
Fax: 978-992-0298

E-mail: info@cellsciences.com
Website: www.cellsciences.com