

## CD274

### Mouse Anti-Human CD274/PD-L1/B7-H1 (22C3) mAb, FITC

<b>Catalog No.</b>	FHJ70121A FHJ70121B	<b>Quantity:</b>	50 tests 100 tests
<b>Alternate Names:</b>	Programmed cell death 1 ligand 1, PD-L1, PDCD1, CD272, B7H1		
<b>Description:</b>	Affinity purified Mouse Anti-Human CD274/PD-L1/B7-H1 (22C3) monoclonal antibody targets Human Programmed cell death 1 ligand 1 (CD274) and was conjugated to FITC.		
<b>UniProt ID (target):</b>	Q9NZQ7		
<b>Target:</b>	Human Programmed cell death 1 ligand 1		
<b>Concentration:</b>	0.5 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Mouse IgG1		
<b>Formulation:</b>	0.01M PBS, pH 7.4, 0.2% BSA, 0.05% Proclin 300		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Clone:</b>	22C3		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Conjugate:</b>	FITC		
<b>Applications:</b>	FCM		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for 12 months. Protect from light. <b>Do not freeze.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

