

## ALK

### Mouse Anti-Human CD246/ALK (ab324) mAb

<b>Catalog No.</b>	FHJ81210A FHJ81210B	<b>Quantity:</b>	50 µg 100 µg
<b>Alternate Names:</b>	ALK tyrosine kinase receptor, Anaplastic lymphoma kinase, CD246		
<b>Description:</b>	Monoclonal Mouse Anti-Human CD246/ALK (ab324) reacts with human ALK tyrosine kinase receptor (CD246)		
<b>UniProt ID (target):</b>	Q9UM73		
<b>Target:</b>	Human ALK tyrosine kinase receptor (CD246)		
<b>Concentration:</b>	1.0 mg/ml, lot specific		
<b>Source:</b>	Mammalian cells		
<b>Isotype:</b>	Mouse IgG1 kappa		
<b>Formulation:</b>	0.01M PBS buffer, pH 7.4		
<b>Purity:</b>	> 95% by SDS-PAGE		
<b>Clone:</b>	ab324		
<b>Purification:</b>	Protein A/G affinity chromatography		
<b>Applications:</b>	FCM		
<b>Storage &amp; Stability:</b>	Store at 2-8°C for up to 1 week, or in working aliquots at -20°C to -80°C for up to 1 year. <b>Avoid freeze/thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.